


face reality™

ACNE CLINIC

How the Acne LED Light Works

LED blue light therapy helps to kill the acne causing bacteria, reduce inflammation and provide a general rejuvenation effect on the skin. Usually, in the case of inflamed acne, topical products will bring down the inflammation rather quickly. However, there are times when the skin needs a little more help in bringing down inflammation. This is where the LED can help.



The combination of blue, red and infrared lights in the LightStim Acne Light get down to the core of the acne causing bacteria. The red lights help to reduce inflammation while at the same time, the blue lights help to kill the acne causing bacteria.

How you will use it

1. Wrap the head of the LED light in plastic (saran) wrap
2. Cleanse your face and dry it
3. Plug in the LED and turn on the switch
4. Place the LED *right on the skin* in the area that has inflammation
5. Leave it on each area for 3 minutes
6. When you have covered the whole area of inflammation, then do the rest of your skincare routine

This is not a substitute for products; nor is it a substitute for your scheduled appointments with us.

****You still need to use your products, see us for regular treatments and get your routine adjusted every two weeks.**

You will need to return this to us within 60 days. We have you sign a credit card receipt for it, but we do not charge you for it. If we do not get it back, we will have to charge you. This LED light is expensive – it is a \$375 investment.

____ I agree to return the Lightstim LED by _____